

Guidelines for succesful installation of IMP software for Geometric Morphometrics

1. Download [mglinstaller.exe](#) and [IMPBasics.exe](#) from the web site at <http://www2.canisius.edu/~sheets/morphsoft.html>

Place these two files in the same folder.

2. Expand (run) [mglinstaller.exe](#) first into a directory.

It will create a number of subdirectories. One of these will be called "bin" and within "bin", there is another subdirectory called "win32".

4. You next need to locate the "win32" subdirectory created when you ran mglinstaller. If you have trouble finding it, use the Windows Find program and search for the file gui_sgl.dll.

5. Expand (unzipp) the compressed IMPBasics.exe into the "win32" directory. The easiest way to run IMP programs is to make sure that all the programs are in the this "win32" folder.

To start any of the IMP programs, go to the "win32" directory and click on the Program icon to start it. That should do it, but if you have problems, don't be afraid to ask.

-Best Wishes

-Dave

H. David Sheets, PhD
Professor
Dept. of Physics
Canisius College
2001 Main St
Buffalo, NY 14208
716-888-2587
sheets@canisius.edu

Adjunct Associate Professor
Dept. of Geology
SUNY at Buffalo
Buffalo NY 14260
(716) 888-2587